

ModuPRO™

Multiple disciplines. One typodont.



At first glance, the ModuPRO may look like your typical typodont that has been around for decades, but it's not!

The ModuPRO has both the upper and lower arches divided into removable segments that house real teeth and artificial teeth. The base fits universally onto many of today's simulation systems.

The ModuPRO is designed to teach multiple dental disciplines in one typodont with added flexibility for the instructor, student and dental board.

Multiple typodonts that are currently used for several disciplines can now be used on one. This increases efficiency and saves the user on the purchase of multiple typodonts.



ModuPRO

Features

Multiple Disciplines

There are currently two types of modules for the ModuPRO – one in which extracted teeth can be placed and one in which artificial teeth can be placed. For example, an entire restorative procedure can be performed with the ModuPRO. Additional disciplines will continue to be added in the future.

Universal fit

The Carrier Trays of the ModuPRO are designed to fit seamlessly in nearly all simulation systems such as pole, magnetic and screw on allowing you to fully utilize your existing simulation lab.

Endo Simulation

Practice canal shaping and access. The ModuPRO Endo allows the student to practice in a simulated patient-based environment. An apex locator can even be used with the ModuPRO Endo Complete.

Easy Radiographs

Each endo module is easy to remove to take multiple radiographic angles, including facial-lingual, mesial distal and a 30 degree angle.

Flexibility

Students can continue to work while the instructor evaluates their work by simply replacing a module rather than the entire typodont.

Examples of ModuPRO Solutions

ModuPRO Endo™ Kit (MPE120 MQR)

The ModuPRO Endo kit comes with all modules to mount extracted teeth, Carrier Trays, Fixing Gel and Apex Putty to mount teeth. This kit is terrific for simulating a variety of endodontic procedures in a patient-based environment.

ModuPRO Pros™ Model (MPR220)

The ModuPRO Pros/Restorative model comes with all modules, Carrier Trays, 28 melamine artificial teeth and adjustable articulator.

** All components also sold separately*

The ModuPRO Provides Savings

Limiting Replacement Costs

	Prosth Dept.	Restorative	Endo Dept.
Entire Typodont Replacement:	Up to \$130	Up to \$130	Up to \$99
Each Module Replacement:	\$16	\$16	\$12.50
Potential Savings \$:	\$114	\$114	\$97

Today, if your typodont becomes damaged, you may have to purchase an entirely new unit, which has little versatility. The ModuPRO, with its modular design, gives you the advantage of only having to purchase a single module/sextant at a time with minimal cost. The long-term cost of replacing typodonts year-to-year will decrease. One module can be replaced, not the entire typodont.

Combining Multiple Disciplines

	Prosth/Restorative/Endo Dept
Today:	Up to \$359 per student
With ModuPRO:	Up to \$229 per student
Potential Savings %:	Up to 36%

** Prices reflect average university pricing*

The ModuPRO has removable sections and is interchangeable between the modules that are used to accommodate extracted teeth as well as the modules for artificial teeth.



To inquire, please contact us at
913-384-7390

5830 Woodson, Suite 3
Mission, KS 66202 U.S.A.

www.acadental.com